



# Carbapenem Susceptibility Testing - Recommendations for Microbiology Laboratories

Utah Department of Health  
Division of Disease Control  
and Prevention  
Bureau of Epidemiology  
HAI Program  
Published August 2014

# Carbapenem Susceptibility Testing - Recommendations for Microbiology Laboratories

---

Prepared by

Felicia Alvarez, MPH  
Healthcare Associated Infections Program; Bureau of Epidemiology  
Utah Department of Health

Special thanks to the following individuals for their subject matter expertise, data resources, editing and consultations.

Primary Children's Hospital  
Judy A. Daly, PhD

Intermountain Healthcare  
Bert K. Lopansri, MD  
E. Kent Korgenski, MS

University of Utah Hospital and Clinics  
Jeanmarie Mayer, MD  
Mark A. Fisher, PhD, D(ABMM)

ARUP Laboratories  
Mark A. Fisher, PhD, D(ABMM)

St. Mark's Hospital  
Joseph Wallis, MT (ASPC)

The Utah Department of Health's *Carbapenem Susceptibility Testing - Recommendations for Microbiology Laboratories* was modeled after the Oregon CRE Toolkit ([http://public.health.oregon.gov/diseasesconditions/communicabledisease/reportingcommunicabledisease/reportingguidelines/documents/cre\\_iguide.pdf](http://public.health.oregon.gov/diseasesconditions/communicabledisease/reportingcommunicabledisease/reportingguidelines/documents/cre_iguide.pdf)) but includes Utah-specific definitions, recommendations and protocols.

Suggested Citation: Utah Department of Health. Carbapenem susceptibility testing - Recommendations for Microbiology Laboratories. Salt Lake City, UT: Utah Department of Health; October 2014.  
[http://health.utah.gov/epi/ARO/CRE\\_Lab.pdf](http://health.utah.gov/epi/ARO/CRE_Lab.pdf)

## Table of Contents

General Recommendations .....	4
Carbapenem-resistant <i>Acinetobacter</i> species (CRAB) Laboratory Response Diagram.....	6
Carbapenem-resistant <i>Escherichia coli</i> and <i>Klebsiella</i> species Laboratory Response Diagram.....	7
Appendix A - List of Utah Laboratories - CRAB/CRE molecular confirmation or Carbapenamase production .....	8
Appendix B - Laboratory Protocol for Detection of Carbapenem-Resistant or Carbapenamase-Producing, <i>Klebsiella</i> spp. and <i>E. coli</i> from Rectal Swabs .....	9
Appendix B – Utah Public Health Laboratory Infectious Disease Test Request Form .....	15

## Recommendations for Microbiology Laboratories

Determine carbapenem susceptibility following the CLSI recommended procedures and interpretive criteria. Between 2010 and 2012, CLSI adjusted susceptibility breakpoints for testing Enterobacteriaceae to carbapenems (see Table below). The 2012 breakpoints increased the sensitivity for carbapenemase detection; laboratories using the 2012 breakpoints do not need to perform a “confirmatory” Modified Hodge Test (MHT) for patient management.

The Utah Department of Health (UDOH) administered a statewide survey in December 2013 and found that most Utah microbiology laboratories used CLSI breakpoints predating the 2010 update and did not perform MHT for confirmatory carbapenemase testing. This is a potential gap in our ability to detect and report carbapenemase-producing CRE (CP-CRE). We recommend that laboratories using pre-2010 breakpoints perform carbapenemase screening and confirmation via MHT. If those laboratories are unable to do their own carbapenemase screening and confirmation, we advise that they send their specimens to a reference laboratory that is able to test (**Appendix A**).

	Breakpoints Predating 2010 Update ( $\mu\text{g/ml}$ ) <sup>1</sup> (through Jan. 2010; M100-S19)			2012 Breakpoints ( $\mu\text{g/ml}$ ) (revised Jun. 2010 and Jan. 2012 <sup>2</sup> ; M100-S22)		
	Susceptible	Intermediate	Resistant	Susceptible	Intermediate	Resistant
Doripenem	n/a	n/a	n/a	$\leq 1$	2	$\geq 4$
Ertapenem	$\leq 2$	4	$\geq 8$	$\leq 0.5$	1	$\geq 2$
Imipenem <sup>3</sup>	$\leq 4$	8	$\geq 16$	$\leq 1$	2	$\geq 4$
Meropenem	$\leq 4$	8	$\geq 16$	$\leq 1$	2	$\geq 4$

In May 2013, the UDOH required the reporting of carbapenem non-susceptible *Acinetobacter*, *E coli*, and *Klebsiella* species. Use the UDOH case definition as follows:

**Confirmed:**

- Isolation of *Acinetobacter* species and:
  - MIC to Imipenem or Meropenem of  $\geq 4 \mu\text{g/mL}$ , or
  - PCR positive for carbapenemase gene
- Isolation of *Klebsiella* species or *Escherichia coli* and:
  - Intermediate or resistant MIC to a carbapenem (refer to table for breakpoints) or
  - Positive Modified Hodge test, or
  - PCR positive for carbapenemase gene

### Screening cultures

CRE screening cultures for outbreak investigations should be performed as recommended by local facility's Infection Prevention and Control staff in consultation with UDOH. The number of surveillance cultures requested will be based on pertinent epidemiology.

- The recommended protocol for screening cultures is included (**Appendix B**). If your laboratory does not have ertapenem (preferred) or meropenem (alternative) disks, contact UDOH HAI epidemiologist. Confirm candidate CRE organisms via routine identification (UPHL should be able to type specimens) and susceptibility (specimens will be sent to a reference lab); keep confirmed CRE isolates until further notification by UDOH.
  - Screening cultures should NOT be billed to the patient; discuss billing with your facility's Infection Prevention and Control department in consultation with UDOH.
  - Facilities may want to consider allocation of 'emergency funds' to pay for active screening cultures in the event of an outbreak or identification of rare organism.
  - Discuss how results of screening cultures will be reported with your facility's Infection Prevention and Control department.

How to send isolates to UPHL for typing only:

Use the General Microbiology Request Form (**Appendix C**).


- In "tests requested", check "other" under isolate identification and write "CRE/CRAB".
- In "comments", please indicate genus and species.
- Send specimen....on a slant; a plate is also acceptable.
- Include collection date, source of specimen, and patient medical record number.

### "Unusual" Isolates


To send any *E. coli*, *Klebsiella spp.*, or *Acinetobacter spp.* or other isolates that are 'unusual\*' to the Utah Public Health Laboratory (UPHL) for further testing at CDC please contact the UDOH HAI Epidemiologist at 801-538-9182. The HAI Epidemiologist will work with you and UPHL to forward the isolate to the Centers for Disease Control and Prevention for further analysis if warranted (not real-time). UPHL will fax results to your laboratory within 3 business days of receipt of test result.

\*Unusual: Rare or unknown isolate.

# Carbapenem-resistant *Acinetobacter* species (CRAB) Laboratory Response Diagram



## Escherichia coli and Klebsiella species Laboratory Response Diagram



**List of Utah Laboratories - CRAB/CRE molecular confirmation or Carbapenamase production**

**LABORATORIES**

Salt Lake City	<p><b>ARUP</b>                      500 Chipeta Way                      Salt Lake City, UT 84108                      801-583-2787  <b>Haleina Muir - Lab Supervisor</b>                      Phone: 583-2787 x 2229</p>
	<p><b>Testing Provided:</b>                      NDM                      KPC</p>
	<p>PCR test</p>
	<p><b>E-Mail:</b> haleina.m.muir@aruplab.com</p>

Salt Lake City	<p><b>INTERMOUNTAIN CENTRAL LABORATORY</b>                      5121 S. Cottonwood St.                      Salt Lake City, UT 84107                      801-507-2244  <b>George Hinde - Lab Supervisor</b>                      Phone: 801-507-2244</p>
	<p><b>Testing Provided:</b>                      NDM                      KPC                      OXA                      IMP                      VIM                      CTXM type ESBL</p>
	<p><b>E-Mail:</b> george.hinde@imail.org</p>

Salt Lake City	<p><b>PRIMARY CHILDRENS HOSPITAL</b>                      100 N. Mario Cappechi Drive                      Salt Lake City, UT 84113                      801-662-2140  <b>Abby Phillips - Lab Supervisor</b>                      Phone: 801-662-2140</p>
	<p><b>Testing Provided:</b>                      KPC</p>
	<p><b>E-Mail:</b> Abby.Phillips@imail.org</p>





## Laboratory Protocol for Detection of Carbapenem-Resistant or Carbapenemase-Producing, *Klebsiella* spp. and *E. coli* from Rectal Swabs

### Purpose

To identify patients colonized with carbapenem-resistant or carbapenemase-producing Enterobacteriaceae in the intestinal tract. Patients who grow these organisms should be placed on Contact Precautions (5) to prevent transmission of the resistant bacteria. The procedure described below is a modification of the procedure described by Landman et al. (4). See the procedural notes for steps in the procedure which can be modified.

### Background

Carbapenem-resistant Enterobacteriaceae (CRE) are usually resistant to all  $\beta$ -lactam agents as well as most other classes of antimicrobial agents. The treatment options for patients infected with CRE are very limited. Healthcare-associated outbreaks of CRE have been reported. Patients colonized with CRE are thought to be a source of transmission in the healthcare setting (1). Identifying patients who are colonized with CRE and placing these patients in isolation precautions may be an important step in preventing transmission (6).

Carbapenem resistance in Enterobacteriaceae occurs when an isolate acquires a carbapenemase or when an isolate produces an extended-spectrum cephalosporinase, such as an AmpC-type  $\beta$ -lactamase, in combination with porin loss. In the United States, the most common mechanism of carbapenem resistance is the *Klebsiella pneumoniae* carbapenemase (KPC).

Detection of carbapenemase production is complicated because some carbapenemase-producing isolates demonstrate elevated but susceptible, carbapenem MICs. CLSI has published guidelines for detection of isolates producing carbapenemases (CLSI document M100) (2). For isolates that test susceptible to a carbapenem but demonstrate reduced susceptibility either by disk diffusion or MIC testing, performing a phenotypic test for carbapenemase activity, the Modified Hodge Test (MHT), is recommended.

Carbapenem resistance and carbapenemase-production in any species of Enterobacteriaceae is an infection control concern. However, the methodology described here focuses on the detection of carbapenem-resistant or carbapenemase-producing *Klebsiella* spp and *E. coli* as this facilitates performance of the test in the microbiology laboratory and, more importantly, because these organisms, especially *Klebsiella* spp. represent the vast majority of CRE encountered in the United States (3).



### **Reagents**

5 ml Trypticase Soy Broth  
10- $\mu$ g carbapenem disks  
MacConkey agar

### **Equipment**

Vortex  
35  $\pm$   $^{\circ}$ C, ambient air

### **Supplies**

100  $\mu$ l calibrated pipettes  
Forceps  
Sterile loops

### **Specimen**

Rectal swab or perianal swab specimen in suitable transport media

### **Special safety precautions**

Biosafety Level 2

### **Quality Control (QC)**

The carbapenem disks that are used in this procedure should be quality control tested using disk diffusion methods and quality control strains as described in the CLSI guideline documents M2 and M100 (2,(2). For example, if the 10- $\mu$ g/mL meropenem disk is used in this procedure, test *E. coli* ATCC 25922 by the disk diffusion method using meropenem disks from the same lot. An acceptable control test will yield a zone size between 28-34 mm. Follow CLSI guidelines for frequency of disk QC testing and corrective action if results are out of range.



## Procedure

Step 1 Day One	<p>Aseptically, place one 10-<math>\mu</math>g ertapenem or meropenem disc in 5 ml trypticase soy broth (TSB) (see procedure note 1)</p> <p>Immediately inoculate the broth with the rectal culture swab</p> <p>Incubate overnight at <math>35 \pm 2^\circ\text{C}</math>, ambient air</p>
Step 2 Day Two	<p>Vortex and subculture 100 <math>\mu</math>l of the incubated broth culture onto a MacConkey agar plate (see procedure note 2)</p> <p>Streak for isolation</p> <p>Incubate overnight at <math>35 \pm 2^\circ\text{C}</math>, ambient air</p>
Step 3 Day Three	<p>Examine the MacConkey agar for lactose-fermenting (pink-red) colonies. More than one colony morphology may represent different species of Enterobacteriaceae (see procedure note 3).</p> <p>It may be necessary to subculture representative colonies of each morphology type to a non-selective media for isolation and/or for susceptibility testing.</p> <p>Screen representative isolated colonies using a phenotypic test for carbapenemase production, such as the Modified Hodge Test (MHT) or test for carbapenem susceptibility using a standardized method and follow the CLSI guidelines for identification of carbapenemase-producing Enterobacteriaceae (see procedure note 4).</p>
Step 4 Day Four	<p>For CRE and/or MHT-positive isolates, perform species-level identification.</p>

## Interpretation/Results

Report all cultures that are positive for CRE or carbapenemase-producing Enterobacteriaceae to the appropriate infection control personnel. Contact Precautions should be implemented for all patients with positive cultures for CRE or carbapenemase-producing Enterobacteriaceae.

## Quality assurance

The ability to recover CRE using this procedure could be assessed as follows: Inoculate 5mL of TSB containing the 10-ug carbapenem disk with a swab that was used to sample a known CRE-negative stool specimen. In addition, inoculate the TSB with 0.5 mL of a  $1 \times 10^5$  CFU/mL suspension of a known carbapenemase-producing isolate (e.g., *K. pneumoniae* ATCC BAA-1705), (see procedural note 5 for suspension preparation) Proceed with Step 2 of the procedure. The carbapenemase-producing *K. pneumoniae* should be recovered on the MacConkey agar.



To test for specificity of the procedure, use a carbapenem susceptible *Klebsiella pneumoniae*, (e.g. ATCC 700603) and follow the same steps. The carbapenem susceptible *K. pneumoniae* isolate should not grow on the MacConkey agar.

### Method limitations

1. Patients may be colonized with CRE or carbapenemase-producing Enterobacteriaceae at a concentration that is not detectable by this method. Studies described by Landman et al. and studies performed at the CDC suggest that the lower limit of detection is between ranges from  $1 \times 10^2$  CFU/ mL to  $1 \times 10^6$  CFU/ mL (4).
2. Non-fermenting gram-negative bacilli with intrinsic mechanisms of carbapenem-resistance, such as *Acinetobacter* spp. and *P. aeruginosa*, will be detected on the MacConkey agar. These isolates should be identified as non-lactose fermenters on the MacConkey agar and therefore would not be picked for characterization. If carbapenem-resistant non-fermenters are present at high concentration, they could overgrow CRE or carbapenemase-producing Enterobacteriaceae on the media and prevent detection of colonization.
3. Enterobacteriaceae can be resistant to carbapenems by mechanisms other than a carbapenemase, the most common of which is expression of an extended-spectrum cephalosporinase, such as an AmpC-type enzyme or an ESBL, combined with porin loss. These isolates will also grow on the MacConkey agar and be identified as carbapenem-intermediate or resistant by standard susceptibility testing but these isolates are negative by the MHT. For isolates that test intermediate or resistant to carbapenems, it may not be necessary to distinguish between these mechanisms of resistance because all carbapenem-nonsusceptible Enterobacteriaceae produce a broad-spectrum  $\beta$ -lactamase, and are therefore an infection control concern. Implementing Contact Precautions for patients colonized with these bacteria would be appropriate. Laboratories may choose to test carbapenem-intermediate or resistant isolates with the MHT to identify carbapenemase-production for epidemiological purposes.

### Procedure notes

1. The procedure described by Landman et al. (4) describes using a 10- $\mu$ g imipenem disk for step 1. However, there are species of Enterobacteriaceae which have intrinsic mechanisms of resistance to imipenem other than a carbapenemase (See CLSI document M100, Appendix G)(2). Therefore, ertapenem or meropenem may provide more specific selection for acquired carbapenem resistance in Enterobacteriaceae.
2. Some laboratories performing cultures for isolation of CRE from rectal specimens place a 10- $\mu$ g carbapenem disk in the first quadrant of the MacConkey plate. Only colonies growing “close” to the carbapenem disk are picked for characterization. No clear criteria for “close” have been established. However, it may be helpful to use either a meropenem or ertapenem disk and then apply the CLSI disk diffusion screening criteria to identify potential carbapenemase-producing isolates (i.e., an ertapenem or meropenem disk zone  $\leq 21$  mm). Note: These zone size criteria



were validated for standardized disk diffusion testing methods as described in CLSI document M2.

3. Carbapenemases are known to exist in several different species of gram-negative bacilli including species of Enterobacteriaceae and *Pseudomonas aeruginosa*. However, carbapenemases are more common in lactose-fermenting species of Enterobacteriaceae (e.g., *K. pneumoniae* and *E. coli*) than in non-lactose fermenting Enterobacteriaceae (e.g. *Serratia marcescens* and some *Enterobacter* spp.) and *P. aeruginosa*. In this procedure, it is suggested that laboratories focus their efforts on detection of resistant lactose-fermenting bacteria to reduce workload. Healthcare facilities that have identified clinical infections with carbapenemase-producing non-lactose fermenting gram-negative species should consider altering this procedure to include characterization of colonies with a morphology that is consistent with those species.
4. The exact procedure for confirmation of CRE or carbapenemase-production should be laboratory-specific and chosen based upon laboratory workflow and the types of isolates causing clinical infections in the patient population served. It may be helpful to refer to the CLSI guidelines for identification of carbapenemase production in isolates that test susceptible to carbapenems in document M100 (2).
5. A  $1 \times 10^4$  CFU/mL suspension of the known carbapenem-resistant or carbapenem-susceptible isolates could be prepared as follows: Dilute 0.1 mL of a 0.5 McFarland standard suspension (equals approximately  $1 \times 10^8$  CFU/ mL), in 9.9 mL sterile water or saline for a 1:100 dilution. From the 1:100, dilute 1.0 mL in 9.0 mL water or saline for a 1:1000 dilution. Add 0.5 mL of the 1:1000 dilution (equals approximately  $1 \times 10^5$  CFU/mL), suspension to the 5 ml TSB for a final concentration of approximately  $1 \times 10^4$  CFU/mL.

## References

1. **Calfee, D., and S. G. Jenkins.** 2008. Use of active surveillance cultures to detect asymptomatic colonization with carbapenem-resistant *Klebsiella pneumoniae* in intensive care unit patients. *Infect. Control Hosp. Epidemiol.* **29**:966-8.
2. **Clinical and Laboratory Standards Institute/NCCLS.** 2009. *Performance Standards for Antimicrobial Susceptibility Testing*. Nineteenth informational supplement. M100-S19. CLSI, Wayne, PA.
3. **Deshpande, L. M., R. N. Jones, T. R. Fritsche, and H. S. Sader.** 2006. Occurrence and characterization of carbapenemase-producing Enterobacteriaceae: report from the SENTRY Antimicrobial Surveillance Program (2000-2004). *Microb. Drug Resist.* **12**:223-230.
4. **Landman, D., J. K. Salvani, S. Bratu, and J. Quale.** 2005. Evaluation of techniques for detection of carbapenem-resistant *Klebsiella pneumoniae* in stool surveillance cultures. *J. Clin. Microbiol.* **43**:5639-5641.



5. **Siegel, J. D., E. Rhinehart, M. Jackson, L. Chiarello, and HICPAC.** 2007. Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings 2007. [www.cdc.gov/ncidod/dhqp/pdf/isolation2007.pdf](http://www.cdc.gov/ncidod/dhqp/pdf/isolation2007.pdf)
6. **Siegel, J. D., E. Rhinehart, M. Jackson, L. Chiarello, and HICPAC.** 2006. Management of Multidrug-Resistant Organisms in Healthcare Settings 2006. [www.cdc.gov/ncidod/dhqp/pdf/ar/mdroGuideline2006.pdf](http://www.cdc.gov/ncidod/dhqp/pdf/ar/mdroGuideline2006.pdf)

# INFECTIOUS DISEASE TEST REQUEST FORM

<b>UTAH PUBLIC HEALTH LABORATORY</b> 4431 SOUTH 2700 WEST TAYLORSVILLE, UTAH 84129 TELEPHONE: (801) 965-2400 FAX: (801) 965-2551 http://health.utah.gov/lab/infectious-diseases	FOR UPHL USE ONLY   LAB# _____ DATE STAMP _____
--	--

PLEASE PRINT CLEARLY FOR ACCURACY.

PATIENT INFORMATION:					
PATIENT STATE OF RESIDENCE:	PATIENT COUNTY OF RESIDENCE:	ZIP CODE:	DATE OF BIRTH (mm/dd/yyyy)	AGE	SEX M F

PATIENT NAME (Last, First):	Is Patient Insured? [ ] Yes [ ] No	STI TESTING ONLY: Is patient MSM? [ ] Yes [ ] No
-----------------------------	---------------------------------------	---

PATIENT ID #	ETHNICITY [ ] Hispanic [ ] Non-Hispanic	RACE [ ] White [ ] Black or African American [ ] American Indian or Alaska Native [ ] Asian [ ] Native Hawaiian or other Pacific Islander
--------------	---	---

<b>PROVIDER INFORMATION</b> Provider Code: _____ Physician: _____ Provider Phone: _____ Provider Email: _____ Secure Fax #: _____	<b>SPECIMEN COLLECTION DATE AND TIME</b> (mm/dd/yy) ____/____/____ Time: _____
--	--

SPECIMEN SOURCE/SITE (CHOOSE 1):			
<input type="checkbox"/> Blood	<input type="checkbox"/> Environmental (specify): _____	<input type="checkbox"/> Plasma	<input type="checkbox"/> Urethra
<input type="checkbox"/> Body Fluid (specify): _____	<input type="checkbox"/> Food (specify): _____	<input type="checkbox"/> Rectum	<input type="checkbox"/> Urine
<input type="checkbox"/> Bronchoalveolar lavage	<input type="checkbox"/> Isolate (source): _____	<input type="checkbox"/> Serum	<input type="checkbox"/> Vagina
<input type="checkbox"/> Bronchial aspirate/wash	<input type="checkbox"/> Lesion (site): _____	<input type="checkbox"/> Sputum (natural / induced)	<input type="checkbox"/> Vomitus
<input type="checkbox"/> Cerebrospinal Fluid	<input type="checkbox"/> Liquid Pap	<input type="checkbox"/> Stool	<input type="checkbox"/> Wound/Abcess
<input type="checkbox"/> Cervix	<input type="checkbox"/> Nasal (aspirate /swab / wash)	<input type="checkbox"/> Throat swab	<input type="checkbox"/> Other (specify): _____
<input type="checkbox"/> (Endo)tracheal aspirate/wash	<input type="checkbox"/> Nasopharyngeal swab	<input type="checkbox"/> Tissue (specify): _____	

BACTERIOLOGY/TUBERCULOSIS TESTS	VIROLOGY / IMMUNOLOGY TESTS
---------------------------------	-----------------------------

<b>Bacteriology Specimen</b> <b>REQUIRED Shipping Temperature:</b> _____ <input type="checkbox"/> Bacterial Culture <input type="checkbox"/> Bacterial ID / Referral Presumptive ID: _____ <input type="checkbox"/> Mycobacterial culture <input type="checkbox"/> Mycobacterial referral Presumptive ID: _____ <input type="checkbox"/> Other (specify): _____	<input type="checkbox"/> C. trachomatis and N. gonorrhoea by NAAT <input type="checkbox"/> Patient is a partner of a 15-24 year old female  <input type="checkbox"/> Herpes/VZV PCR (HSV-1, HSV-2, VZV)  <input type="checkbox"/> Virus Identification (culture) Virus suspected _____  <input type="checkbox"/> Cytomegalovirus  <input type="checkbox"/> Varicella zoster virus  <input type="checkbox"/> Multi-Pathogen Respiratory Panel (Includes Adenovirus, Coronavirus, Human Metapneumovirus, Rhino/Enterovirus, Influenza A, Influenza B, Parainfluenza 1-4, RSV, Bordetella pertussis, C. pneumoniae, M. pneumoniae)  <input type="checkbox"/> Influenza A & B virus PCR (with subtyping) <input type="checkbox"/> Hospitalized w/ Influenza-like illness <input type="checkbox"/> Other (i.e., cluster investigation) Cluster location: _____ Other reason for testing: _____  <input type="checkbox"/> West Nile virus IgM (Human)	<input type="checkbox"/> QuantiFERON-TB Gold <b>REQUIRED information:</b> Blood draw date/time: _____ Incubation at 37°C completed? [ ] Yes [ ] No Signature: _____ Incubation start date/time: _____ Incubation end date/time: _____  <input type="checkbox"/> Syphilis IgG EIA (includes confirmatory testing) <input type="checkbox"/> RPR (suspect acute infection/previous positive)  <input type="checkbox"/> HIV Antigen/Antibody (includes confirm. testing) <input type="checkbox"/> Previous positive  <input type="checkbox"/> Hepatitis C Antibody <input type="checkbox"/> Add HCV RNA Testing if Positive  <input type="checkbox"/> Hepatitis C RNA (Qualitative; Antibody screen not included)  <input type="checkbox"/> Hepatitis B Antibody  <input type="checkbox"/> Hepatitis B Antigen  <input type="checkbox"/> Hantavirus (Sin Nombre) IgG/IgM <input type="checkbox"/> Acute Serum (mm/dd/yy) ____/____/____ <input type="checkbox"/> Convalescent serum (mm/dd/yy) ____/____/____
<b>BIOTERRORISM TESTS</b> <u>(Notify Lab before submitting)</u> <input type="checkbox"/> Bacillus anthracis Detection/Identification <input type="checkbox"/> Brucella species Detection/Identification <input type="checkbox"/> Brucella antibody <input type="checkbox"/> Burkholderia mallei/pseudomallei Detection/ID <input type="checkbox"/> Clostridium botulinum culture & toxin <input type="checkbox"/> Coxiella burnetii Detection <input type="checkbox"/> Francisella tularensis Detection/Identification <input type="checkbox"/> F. tularensis antibody <input type="checkbox"/> Orthopox viruses Detection Virus Suspected: <input type="checkbox"/> Vaccinia virus <input type="checkbox"/> Varicella zoster virus <input type="checkbox"/> Variola virus <input type="checkbox"/> Yersinia pestis Detection/Identification <input type="checkbox"/> Yersinia pestis antibody <input type="checkbox"/> Other (specify): _____		

ADDITIONAL INFORMATION	
[ ] Other Disease Suspected: _____	[ ] Referral Test to CDC (form <b>REQUIRED</b> ) specify: _____ Contact UPHL for CDC form

COMMENTS:

**Utah Department of Health  
Healthcare-Associated Infections  
Prevention & Reporting Program**  
P.O. Box 142104  
Salt Lake City, UT 84114  
<http://health.utah.gov/epi/diseases/HAI/index.html>  
PH 801-538-6191 • FAX 801-538-9923

